

## Development and analytical validation of a fully-automated platform for quantification of MetaSites to predict systemic metastasis.



Mike Peterson<sup>1</sup>, Lonnie Graves<sup>1</sup>, Michael J. Donovan<sup>2</sup>, Douglas Hamilton<sup>1</sup>, Olga Pozdnyakova<sup>3</sup>, and Mark Gustavson<sup>1</sup>  
<sup>1</sup>MetaStat, Inc. (Boston, MA), <sup>2</sup>Mt. Sinai Medical Center (New York, NY), Brigham and Women's Hospital (Boston, MA)

## Abstract

**Background:** Tumor metastasis is responsible for the majority of solid tumor related deaths. Diagnostic accuracy that is accurate and predictable risk of cancer dissemination can provide useful information for optimizing personalized treatment strategies. A previous study demonstrated the use of immunohistochemical methods, the number and distribution of Metastases in the human microenvironment has significant predictive value for patients with B-positivity, early stage invasive breast cancer. A Metastis-e is the x-ray projection of a multicellular structure comprised of a blood vessel, the macrophage immune cell, and an invasive tumor cell. To reduce pathologist variability, an objective and reproducible automated laboratory process workflow to identify and quantify Metastases in a clinical gray area was developed. **Methods:** Digital pathology image coupling with image analysis tools was employed to develop a fully automated, objective method and workflow for quantification of Metastases in formalin fixed paraffin embedded tumor samples. Using this method, areas for analysis and quantification of Metastases were automated by integrating high resolution automated microscopy with multiple image analysis algorithms. A pathologist ensured overall diagnostic quality of the samples in addition to approving individual images for Metastase scoring. **Results:** In this analytical validation study, the platform was demonstrated to be greater than 9% reproducible with a mean coefficient of variation of 6.6% ( $n=35$ ) for 3 independent measurements of the same slide. Further, Metastase scores showed correlation coefficients (Pearson's R) greater than 0.98 between measurements with no significant difference in absolute values by repeated measures analysis. Importantly, Metastase scoring on independently stained tumor sections showed greater than 90% reproducibility, indicating minimal heterogeneity within the tumor with respect to Metastase score, section to section. Additionally, daily day-to-day Metastase scores showed correlation coefficients (Pearson's R) greater than 0.90 between staining sections with no significant difference in mean Metastase scores. **Conclusion:** Taken together, these data demonstrate the successful development and analytical validation of a fully automated, highly reproducible Metastase quantification platform. With further development and analytical validation of this test, it is now possible to provide physicians with information regarding the aggressiveness of patient tumors and accurate prediction of cancer metastasis. This method is being further validated in a large ( $n=481$ ) case control clinical study.

## Introduction

Metastatic breast cancer is a driver-based triple-immunotherapy to hematologic tumor cells that identify the site of intratumoral angiogenesis and predict the well-elliptical metastatic tumor cells in early stage non-regressive and non-progressive invasive breast cancer. Breast cancer is one of the most frequently diagnosed cancers in women with an estimated 20,800 new cases diagnosed in the US in 2018. Aggressive surgery leading to metastasis is responsible for the majority of the estimated 40,370 deaths in 2018. A combination of immunotherapy with chemotherapy is the standard of care for metastatic breast cancer.

30,760 deaths (95% CI 28,222-33,308). Approximately 30,000 women diagnosed with early stage breast cancer who undergo primary treatment for their disease will die from metastatic disease over the course of their lifetime. Although the overall prognosis is a favorable one for all patients with early stage, estrogen receptor-positive invasive breast cancer following primary therapy, women with advanced disease have a poor prognosis. Women with metastatic disease have a median survival of approximately 2 years, and the 5-year survival rate is approximately 15%. The prognosis for women with metastatic disease is significantly worse than for women with early stage disease, and the overall survival rate is approximately 10% at 5 years. The prognosis for women with metastatic disease is significantly worse than for women with early stage disease, and the overall survival rate is approximately 10% at 5 years.

In current issue of *Mobile Health*, the WeSite module quantified reproducibly by a protocol called object counting. This was been enabled by the fast counting of cellular structure counts per unit area and interquartile variability ( $2.7 \pm 0.6\%$  coefficient of variation) in relation to counting method, specificity for quantification in the absence of ablation. Mendeisite module resulted in only 75% interobserver reproducibility. In order for WeSite module to be deployed in the clinical laboratory setting and to meet these standard required for clinical use it will be necessary to develop reproducible process for quantification of WeSite instrument. Hence we demonstrate the development and validation of a robust and objective method for WeSite quantification using image analysis and bioinformatics.

acquired (up to 100) frames of  $\lambda$  images to run using custom software in Form (Pernix Image) may also align them (See Figure 1c). Individual metastases were identified using Form and Visopharm (Visopharm, Hornsholm, Denmark) image analysis software. In Form, areas of at least 10% of the total image representing chromogen-specific channels for CD31, C068, and Pan Melo. These images were analyzed using the Visopharm Metastate identification algorithm where individual metastases were identified as structures which conform to established criteria.

**Pathology QC**  
All samples in a study were reviewed by a pathologist and any images not meeting standard for invasive breast cancer were removed from analysis. Furthermore, the pathologist will review each slide for quality and specific staining patterns. Two pathologists independent of QCP processes checked 97.30% of cases in our study.


cases, final n=8). Analytical precision was defined as repeated analysis ( $n=8$ ) of the same slide through (3) times with a primary a priori determined bench mark of 95% reproducibility by ROC analysis. Analytical performance was defined as completely independent staining and analysis ( $n=5$ ) on serial tissue sections on a bone specimen at a day with a primary a priori determined bench mark of 90% reproducibility by ROC analysis. Analytical accuracy was defined as comparison of automated MIF-100 quantification with the academic assay.

determined benchmark of 98.0% reproducibility by ROC analysis. Analytical accuracy was defined as comparison of automated MetStat quantification with the academic assay reference standard with a primary a priori determined benchmark of 98.0% accuracy by ROC analysis. All statistical analysis was performed in SPSS v. 22 (SPSS Inc., Chicago, IL).

- 1.) <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3907160/>
- 2.) Toss et al. Molecular characteristics and targeted therapeutic approaches in breast cancer. *Breast Cancer Res.* 2015 Apr 23;17(6).
- 3.) Rohan et al. Tumor Microenvironment of Metastasis and Risk of Distant Metastasis of Breast Cancer. *JNCI J Natl Cancer Inst.* 2014 Jun;106(7)


4.) Polley et al. An international K6 reproducibility study. J Natl Cancer Inst. 2013 Dec 18;105(24):1897-906

## Results




**Figure 1:** Auto mated field capture and image analysis for objective quantification of Me<sub>3</sub>Ts in breast cancer samples. (A) Inform coupled with the Vectra 2 microscopy system is used to automatically take high resolution 20X fields of view (FOVs) to identify Me<sub>3</sub>Ts. (B) Example of high resolution 20X FOV showing Me<sub>3</sub>Ts positive tumor (red), CD31-positive blood vessels (blue) and CD68-positive macrophages (brown). (C) Visiopharm image analysis software is used to identify tumor (red), vessels (blue), and macrophages (brown) and the n (D) Me<sub>3</sub>Ts are identified by determining percentage of juxtaposed cellular interfaces (blue/red interface).

Analytical Precision and Performance assessment was performed by comparing reproducibility of MetaSite™ Breast scores for multiple sites (3) independent samples/analysis of the same slide (Precision) or three (3) independent staining runs and analysis on separate tissue sections (Performance). Figure 2 shows linear regression and mean analysis for continuous MetaSite™ Breast score data. Strong positive Pearson's ( $R = 0.99$ ) and slope approach 1.00 ( $0.99 \pm 0.01$ ) were observed in all comparisons with no statistically significant difference in continuous MetaSite™ Breast scores across all treatments (precision) or staining runs (Performance). Figure 3 shows ROC analysis categorical (Low vs. Intermediate/High) Low/Intermediate v. High (High KU point) MetaSite™ score comparisons. Area Under the Curve (AUC) measurement is indicating degree of agreement ranged from 0.57-0.99 for precision and 0.93-0.96 for performance with an overall Pearson Correlation Coefficient of 0.66.




**Figure 2:** Analytical Precision and Performance Validation. Linear regression and means analysis for (A) Analytical Precision (repeat analysis is the same slide) and (B) Analytical Performance (staining and analysis of serial sections on separate days). Shown are XY Scatterplots (Repeat/Day 1 v. Repeat/Day 3 only) and Means Plot (inset) showing very strong correlation (Pearson's R) and slopes approaching 1 with no significant difference in mean scores by repeated measures analysis. Taken together these data demonstrate a very high level of precision and performance for raw continuous MetaSite scores.

## Results Continued



**Figure 3:** Analytical Precision and Performance Validation continued. ROC analysis for (A) Analytical Precision (repeat analysis of the same slide) and (B) Analytical Performance (staining and analysis of serial sections on separate days) using low and high [inset] cutpoints set by taking tertiles and comparing to subsequent repeat analyses. Shown are comparisons for Repeat Day 1 versus Repeat Day 3 with indicated area under the curve (AUC) measurements. Taken together these data demonstrate a high level of reproducibility with respect to group classification both for repeated analysis (Precision) and different sections of the same tumor sample (Performance).

Analytical Accuracy Assessment was performed by comparing automated MetaSite™ Breast scores with the reference standard, semi-manual pathologist TME counting on the same slide. ROC analysis (Figure 4) shows a high level of concordance, 84% for the low cutpoint (low v. intermediate/high) and 90% for the high cutpoint (low/intermediate v. high).



**Figure 4:** Accuracy Analysis Validation. Continuous vs automated MetaSite scoring was compared to semi-manual *pa*thologist scoring (reference standard) by ROC analysis using the original Pathologist cutpoints (Rohan *et al.*)<sup>1</sup> for (A) pathologists low v. intermediate/high risk cutpoint and (B) pathologists low/intermediate v. high risk cutpoint *ts*. Data demonstrate a high degree of concordance showing 84% and 90% concordance respectively.

## Conclusio

- The MetaSite™ Breast assay, utilizing automated digital pathology methods, demonstrates high reproducibility, accuracy and analytical precision, thus enabling MetaSite™ Breast to be used in the clinical setting.
  - With precision and performance of the assay, specifically %CV < 10%, the MetaSite™ Breast test meets or exceeds current industry standards for analytical performance of *in situ* tissue-based diagnostic tests.

Contact Mark Gustavson, PhD for more information or reprints ([mgustavson@metastat.com](mailto:mgustavson@metastat.com))