



Aurora Masks by Inogen®

Aurora N1 Nasal Mask

The Aurora N1 Nasal mask offers patients a secure fit in a minimalist cushion design without the compromise of comfort or performance.

In a 90-day study of CPAP nasal cushion mask users (n=44), a majority reported being highly satisfied with the Aurora N1 mask.¹



Aurora N1 Nasal Mask



In the 90-day study (N=44 nasal cushion users), 81% found the Aurora N1 mask to be just as or more durable as their current mask.¹

1 High-Quality Ergonomic Crown-Shape Headgear
Mask headgear design features a triple-material design with ultrasonic edge sealing technology, promoting comfort and durability.

2 Adaptive Silicone Cushion
Medical-grade silicone cushion design features an adaptive design with a range of up to 15 degrees at the nose bridge, adjusting to different face shapes.

3 Quiet Honeycomb Vent Design
Honeycomb vent design reduces airflow intensity to provide gentle and quiet exhalation for patients, designed to improve comfort and adherence.

4 Quick-Release Headgear Clips
Quick-release headgear clips allow for ease of putting on and taking off the mask, making adjustments easy and user-friendly.

SPECIFICATIONS

Therapy Pressure Range	4-30 cmH ₂ O
Tube Connection	22mm
Flow Resistance (At 50 L/min)	≤0.5 cmH ₂ O
Flow Resistance (At 100 L/min)	≤1.5 cmH ₂ O
Cushion Sizes	S/M/L/XL

provider.inogen.com

1. AMG CPAP Mask Trial Study Aug 2025, 90-day study of 44 patients using CPAP therapy with a nasal cushion mask for: at least 6 months, at least 5 nights a week, and at least 4 hours a night.

Inogen® Aurora Masks

The mask is a noninvasive accessory used for channeling airflow (With or without supplemental oxygen) to a patient from a positive airway pressure (PAP) device such as a continuous positive airway pressure (CPAP) or bi-level system. For safe and proper use of the Inogen® Aurora masks (including intended use, contraindications, and warnings), please refer to the instructions for use at www.inogen.com/auroramasks.

PM250183 EN_EX_USA | Rx Only. © 2026 Inogen, Inc. | 859 Ward Drive, Suite 200, Goleta, CA 93111

Inogen® and other related product names contained herein are trademarks of Inogen, Inc. The usage of any Inogen, Inc. trademark is strictly forbidden without the prior consent of Inogen, Inc. All other trademarks are trademarks of their respective owners or holders.