



Aurora Masks by Inogen®

Aurora P1 Nasal Pillows Mask

The Aurora P1 nasal pillows mask is easy for patients to use while delivering performance in a compact design.

In a 90-day study (N=43 nasal pillow mask users), a majority of patients reported the Aurora P1 mask is just as or more comfortable than their current mask.¹



Aurora P1 Nasal Pillows Mask



In the 90-day study (N=43 nasal pillow users), 95% found the Aurora P1 mask to be just as or more durable as their current mask.¹

- 1 Durable Headgear Design**
Mask headgear design features a triple-material design with ultrasonic edge sealing technology, promoting wearing comfort and product lifespan.
- 2 Comfortable Silicone Cushion**
Medical-grade rebound silicone cushion is soft, airtight, and hypoallergenic-designed to reduce facial pressure and improve patient fit and comfort.
- 3 Quiet Honeycomb Vent Design**
Honeycomb vent design reduces airflow intensity to provide gentle and quiet exhalation for patients, designed to improve comfort and adherence.
- 4 360° Swivel Elbow**
360° Swivel elbow allows for flexible movement during sleep, adapting to different sleeping positions effortlessly.

SPECIFICATIONS

Therapy Pressure Range	4-20 cmH ₂ O
Tube Connection	22mm
Flow Resistance (At 50 L/min)	≤1.0 cmH ₂ O
Flow Resistance (At 100 L/min)	≤3.5 cmH ₂ O
Cushion Sizes	S/M/L

provider.inogen.com

1. AMG CPAP Mask Trial Study Aug 2025, 90-day study of 43 patients using CPAP therapy with a nasal pillow mask for: at least 6 months, at least 5 nights a week, and at least 4 hours a night.

Inogen® Aurora Masks

The mask is a noninvasive accessory used for channeling airflow (With or without supplemental oxygen) to a patient from a positive airway pressure (PAP) device such as a continuous positive airway pressure (CPAP) or bi-level system. For safe and proper use of the Inogen® Aurora masks (including intended use, contraindications, and warnings), please refer to the instructions for use at www.inogen.com/auroramasks.

PM250184 EN_EX_USA | Rx Only. © 2026 Inogen, Inc. | 859 Ward Drive, Suite 200, Goleta, CA 93111

Inogen® and other related product names contained herein are trademarks of Inogen, Inc. The usage of any Inogen, Inc. trademark is strictly forbidden without the prior consent of Inogen, Inc. All other trademarks are trademarks of their respective owners or holders.